



Nursing Bottle Brush for Cordless Drill



Shown with cordless power drill
(not included)

Used with a rechargeable power drill, this brush makes quick work of cleaning large numbers of two, three and four quart nursing bottles. Brushes will reach every corner and bend of the bottle. We recommend running the drill at low to medium RPM for maximum life of the brush. 18-1/2" overall length, 5-1/2" wide. Wt. 6 oz.

- 1030-147** Soft Bristle Drill-Operated Nursing Bottle Brush.
- 1030-148** Stiff Bristle Drill-Operated Nursing Bottle Brush.



The Coburn Company, Inc. • P O Box 147 • Whitewater, WI 53190 U.S.A.

800 776 7042

262 473 2822


coburn.com

info@coburn.com

E-Z Nurse™ Plastic Nursing Bottles



996
Three Quart
Calf Nurser



Nursers are available with standard screw-cap and nipple or with one-piece snap-on nipple. The snap-on nipple will also fit the Carnation bottle. Mouth of bottles measure 2-5/8" OD.

Two-Quart Nursers

- 911 Complete Nurser including screw cap, calf nipple and two-quart bottle. Carton of 12 weighs 6.5 lbs.
- 1030-203 Complete Nurser including snap-on calf nipple and two-quart bottle.
- 913B Repl. 2-quart bottle only, screw cap type.
- 919B Repl. 2-quart bottle only, snap-on type.
- 915R Wire Bottle Rack for two-quart nursers. Carton of 12 weighs 7 lbs.

Three-Quart Nursers

- 996 Complete Nurser including screw cap, calf nipple and three-quart bottle. Ctn. of 6 weighs 4.6 lbs.
 - 1030-204 Complete Nurser including snap-on calf nipple and three-quart bottle.
 - 918R Wire Bottle Rack for three-quart nursers. Carton of 12 weighs 9 lbs.
 - 999 Replacement 3-quart bottle, screw cap type.
 - 918B Replacement 3-quart bottle, snap-on type.
 - 846 24" Bottle Cleaning Brush for 3-quart bottles.
 - SP1515 24" White Bottle Brush for 3-quart bottles.
- ### Replacement Nipples & Caps
- 914N Replacement screw-on type calf nipple.
 - 916C Replacement screw cap.
 - 1030-201 Replacement snap-on type calf nipple.

E-Z Nurse™ "Maximum" Calf Nursers



- "Maximum" Snap-On Calf Nipples provide maximum bottle grip and maximum nipple length. They are perfect for hutches. Nipples provide increased milk flow and are "E-Z" to clean.
- 994 MAX 2 Qt. Nurser with Snap-On Nipple. Wt. .4 lb.
 - 995 MAX 3 Qt. Nurser with Snap-On Nipple. Wt. .6 lb.
 - 993 Repl. MAX Snap-On Nipple. Wt. 1.4 oz.



Nursing Bottles




- Nursing Bottles are made of strong, durable polyethylene plastic. Snap-on or screw-on nipples are pliable and long lasting. Screw-on nipple has vacuum air vent. Mouth of bottle measures 2-1/4" OD.
- 9812 2 Qt. Nurser with Screw Cap & Nipple.
 - 9312 2 Qt. Nurser with Snap-On Nipple.
 - 93CB 2 Qt. Bottle only for screw-on nipples.
 - 98CB 2 Qt. Bottle only for snap-on nipples.
 - 9313 3 Qt. Nurser with Screw Cap & Nipple.
 - 9813 3 Qt. Nurser with Snap-On Nipple.
 - 9313B 3 Qt. Bottle only for screw-on nipples.
 - 9813B 3 Qt. Bottle only for snap-on nipples.
 - 935CN Repl. Screw-on Calf Nipple.
 - 935C Repl. Screw-on Cap.
 - 98CN Repl. Snap-On Calf Nipple.



POLY DOME Jumbo 3-Qt. Nursing Bottles

- More and more veterinarians and universities say that 2 quarts of milk is not enough in cold weather. In fact, some vets claim it is not enough under any conditions. For faster gain and healthier calves, switch to Poly Dome 3 quart nursing bottles. Easy to wash and sterilize. Mouth of bottle measures 2-1/4" OD. Rack and carrier, available separately, are made of powder-coated 3/16" steel wire and reinforced for strength.
- PD3 3 Quart Bottle with snap-on nipple.
 - PD4 Replacement Nipple.
 - PD5 Replacement 3 Quart Bottle.

Manna Pro 2 Qt. Nursing Bottles



- Manna Pro calf nursing bottle is made of durable plastic for long-lasting performance and ease of cleaning and sanitizing. The bottle is designed for "heads up" feeding, the natural way. Bottle holds 2 quarts of milk replacer. Includes bottle and snap-on nipple.
- 0942029906 Manna Pro Nurser for Calves.
 - 0942059901 Repl. Calf Nipple.



Nursing Bottle



- Milk Bar Nursing Bottle features the exclusive Milk Bar slow-nursing nipple. Saliva contains enzymes for digestion and antibiotic properties for boosting calves' immunities and slow-nursing helps maximize saliva production. Calves have a strong sucking instinct and other nipples may provide too much milk too fast which can choke the calf and cause milk to enter the lungs or the wrong stomach where milk scours can result.
- 265-9001 Two-Quart Nursing Bottle with Milk Bar Snap-On Nipple.
 - 265-9002 Two-Quart Nursing Bottle with Milk Bar Snap-On Nipple for colostrum feeding, yellow.



Snap-On Nipples for Nursing Bottles



- Designed to replicate nature, the Milk Bar Nipple provides milk at the proper rate to stimulate saliva production which is essential for nutrition and for the development of the calf's digestive system. Milk Bar Nipples fit our snap-on style nursing bottles or any bottle with top diameter of 2-1/2" to 2-5/8" (57 to 67cm). Nipples are available bulk or package of 2. Yellow nipple for colostrum is made of a softer rubber formulation.
- 265-0005 Milk Bar Snap-On Nipple, black, bulk.
 - 265-0005-2 Milk Bar Snap-On Nipple, black, pkg. of 2.
 - 265-0006 Milk Bar Snap-On Nipple for colostrum feeding, yellow, bulk.



The Coburn Company, Inc. • P O Box 147 • Whitewater, WI 53190 U.S.A.

800 776 7042

262 473 2822

coburn.com

info@coburn.com